

# Cuna

Gently supports forward, upright, and backward leaning postures. A chair that sways comfortably.



Inorganic Antibacterial Agent  
JP0122675X00041

\*Inhibit proliferation of bacteria on product surface.  
- Backrest

# cuna

## ■ Selection

K13-ZZ11C-BKGHM61

Black body (BK)  
Upholstered seat – Dark greige (GHM6)



K13-ZZ11C-E2GH1K1

Light grey body (E2)  
Upholstered seat – Light greige (GH1K)



K13-ZZ11C-E2GH2J1

Light gray body (E2)  
Upholstered seat – Opal green (GH2J)



K13-ZZ11C-E2GH1S1

Light gray body (E2)  
Upholstered seat – Olive yellow (GH1S)



※ Prices if shown: List prices (Excluding Tax)

## Supports optimal posture by natural swaying

The inclined leg frame transfers body weight while swaying to prompt comfortable seating. Supports forward, upright, and backward leaning postures to aid the working posture. A chair for the body to sway and move naturally, reducing pressure strains on the body from long-hour use.

Backward leaning



Forward leaning




Forward leaning posture suited for focus work



Refresh while swaying

## ■ Seat cushion ideal for long-hour use + Contoured arms for ease-of-use



### Seat cushion structure

Plano-convex shaped seat with firmer middle core. The molded urethane foam cushion completes seating comfort suitable for everyday high-traffic use.



### Antimicrobial agent kneaded PP backrest\*

The growth of bacteria is prevented by 99% or more for 24 hours (in comparison to non-antibacterial treated surfaces). Kohkin (inhibiting bacteria) as tested in accordance with ISO22196 International standard.

\* Antimicrobial performance is for the backrest only. Antimicrobial treated area reduces specific bacteria on the surface and not a means to guarantee the prevention of infection. Further, the antimicrobial treatment is not intended to treat or prevent infectious diseases.

## Color Codes

Body (Poly)	BK	E2	E2	E2
	Black	Light gray	Light gray	Light gray
Seat (Fabric)	GHM6	GH1K	GH2J	GH1S
	Dark greige	Soft beige	Opal green	Olive yellow
Color code (for use in product code)	BKGHM6	E2GH1K	E2GH2J	E2GH1S